made in GERMANY

ORIHOMEDICAL
G m b H · I M P L A N T S

MODULAR KERRISON HYBRID SYSTEM

NOTIZEN / NOTES

INHALTSVERZEICHNIS / INDEX

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MODULAR KERRISON HYBRID SYSTEM

Since the first publication of Percutaneous operations of Lumbar intervertebral discs in the early 1970s, the full-endoscopic surgical technique has been developed rapidly. Since that time, a wide range of new surgical and access techniques have been developed.

Due to the increasing miniaturization in endoscopy, there is a growing demand for delicate but even robust and easy to clean hand instruments. All these requirements are completely fulfilled by our new modular Endo Kerrison Hybrid System.

This unique hybrid product line combines the characteristics:

- Modularity and Robustness
- Disassembly- and Cleanability
- Ergonomics and Individuality

better than any other system.


Developed and made in Germany, available in various shaft dimensions, shaft designs and working lengths, the ENDO-Kerrison punches are compatible with all established endoscopic systems on the market.

The individual combinability with the two types of handles made of medical-grade stainless steel, and optional with silicone coating, guarantees optimal handling for every user.


Thanks to the unique design, which allows the complete disassembling of the Endo-Kerrison punch, all requirements for cleaning in accordance with EN 17664 are met and the safety of patients and users is guaranteed.

Excellent cutting performance thanks to the large cutting range

Perfect to clean integrated flushing channel


Optimum cutting geometry 90° and 40° (130°) Cutting Edge

*Take-apart design...
allows the alternation
between various shaft designs*

*Open construction allows
optimal cleanability*


*Perfect ergonomics due
to the adjustability of
the cutting edge. Can be
rotated around 360° and
adjusted in four positions.*

*Ergonomics and High Grip!
Ergonomic handle made
of medical-grade
stainless steel, if desired
with silicone coating.*

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Endo - Kerrison Modular Handle
with silicone coating blue (Cat. No.55-000.00)

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Endo - Kerrison Shaft System
Variants overview at page 10 - 11

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Kerrison - Shaft System (short)
Variants overview at page 13 - 14

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Endo - Kerrison - lateral relieved
Variants overview at page 12

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INSTRUCTION FOR USE AND ASSEMBLING FOR MODULAR KERRISON HYBRID SYSTEM

Intended Use

The punches are designed for adequate use only. The decision on the suitability of the punches for the intended application lies with the user. We assume no liability for the consequences of improper application, maintenance or preparation of lubrication cleaning, reprocessing and sterilization.

Warning / Caution Notice

Twisting and levering of the punches must be avoided, long shafts – punches with 290 mm and 350 mm working length are the intended exclusive use for minimal invasive endoscopic application – indication, excepted are only standard shafts – punches with working length 220 mm. Mechanical overloading can damage the instrument or even cause break or deformation of the jaw.

We decline any liability for damages, evidenced to mechanical overloading!

Warranty

We guarantee that our products are free from any manufacturing processing defects.

Functional Testing

New instruments must undergo a thorough visual and functional inspection before us.

Repairs

Even under appropriate use of surgical instruments are subject to more or less severe wear.

Servicing and repair may only be carried out by our self as the manufacturer. Surgical instruments which must be send back for repair have to be cleaned, disinfected in accordance to reprocessing requirements of EN17664, also must be returned in sterilized condition.

Preparation Reprocessing Instructions

Please observe our separate cleaning -and reprocessing instruction and in accordance to reprocessing requirements of EN17664. All thorough cleaning of all parts is indispensable. Blood, pus and protein residues can cause damages of the material and function.

The punch shafts are easily be disassembled for cleaning and disinfection, the modular punch handle which must be not disassembled, only remove the blue shaft connecting nut before reprocessing, qualifies and applies a save cleaning and disinfection too.

After cleaning and disinfection of all moveable surgical instrument parts have be lubricated with a licenced lubrication oil for surgical instruments only.

Sterilization

All instruments - components are supplied non sterile and must be cleaned sterilized before use.


Intended Purpose

Surgical instruments for micro surgical decompression of lumbar spinal stenosis.

GEBRAUCHS -UND MONTAGANLEITUNG FÜR MODULARES KERRISON HYBRID SYSTEM | INSTRUCTION FOR USE AND ASSEMBLING FOR MODULAR KERRISON HYBRID SYSTEM

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

ORIHOMEDICAL


G m b H · I M P L A N T S

THE SPINE ENDOSCOPE THE PERFECT ENDOSCOPIC VIEW

SOLUTIONS FOR SPINE SURGERY

THE SPINE ENDOSCOPE


Features

- Sapphire glass window
- Gold soldered opto-mechanics
- Large working channels
- Incorporated irrigation and suction channels
- Stainless Steel
- Fully autoclavable

Specification/ Item Number	30-0734-00	30-0780-00
Viewing Direction	30°	30°
Outer Diameter	Ø6,3 mm	Ø6,9 mm
Working length	193mm	193mm
Total length	260mm	260mm
Working channel	Ø3,75 mm	Ø4,1 mm
Suction/ irrigation channel	2x Ø1,4mm	2x Ø1,4mm

Artikel-Nr. Cat. No.	Bezeichnung Description
30-0734-00	Spine Endoscope Ø 6.3 mm L= 193mm, 30° Autoclavable 134°
30-0780-00	Spine Endoscope Ø 6.9 mm L= 193mm, 30° Autoclavable 134°

LITERATURQUELLEN - BIBLIOGRAFIE | LITERATURE SOURCE - BIBLIOGRAPHY

1. Arnoldi CC, Brodsky AE, Cauchoix J et al. Lumbar spinal stenosis and nerve root entrapment syndroms: Definition and classification. *Clinical Orthopaedics* 1976; 115: 4–5
2. Boden S, Davies DO, Dina TS et al. Abnormal magnetic resonance scans of the lumbar spine in asymptomatic subjects. A positive investigation. *Journal of Bone Joint Surgery American* 1990; 72: 403–408
3. Caspar W. Technique of microsurgery. In: Williams RW, McCulloch JA, Young PH (eds) *Microsurgery of the lumbar spine*. Aspen: Rockville, 1990: 105–122
4. Deen HG Jr, Zimmermann BS, Swanson SK, et al. Assessment of bladder function after lumbar decompressive laminectomy for spinal stenosis: A prospective study. *Journal of Neurosurgery* 1994; 80: 971–974
5. Deyo RA, Cherkin DC, Loeser JD et al. Morbidity and mortality in association with operations on the lumbar spine: The influence of age, diagnosis, and procedure. *Journal of Bone Joint Surgery American* 1992; 74: 536–543
6. Epstein JA, Epstein BS, Rosenthal AD, et al. Sciatica caused by nerve root entrapment in the lateral recess: The superior facet syndrome. *Journal of Neurosurgery* 1972; 36: 584–589
7. Epstein NE. An evaluation of varied surgical approaches employed in the management of 170 far lateral lumbar disc herniations. *Journal of Neurosurgery* 1995; 83: 648–656
8. Epstein NE. Symptomatic lumbar spinal stenosis. *Surgery of Neurology* 1998; 50: 3–10
9. Epstein NE, Epstein JA. Lumbar decompression for spinal stenosis: Surgical indications and techniques with and without fusion. In Frymoyer JW (ed). *The adult spine: Principles and practice, volume 2*, 2nd. Edition New York: Lippincott-Raven, 1997: 2055–2088
10. Epstein NE, Epstein JA, Carras R, et al. Coexisting cervical and lumbar spinal stenosis: Diagnosis and management. *Neurosurgery* 1984; 15: 489–496
11. Hart LG, Deyo RA, Cherkin DC. Physician office visits for low back pain. Frequency, clinical evaluation, and treatment patterns from a U.S. national survey. *Spine* 1995; 20: 11–19
12. Katz JN, Lipson SJ, Lew RA, et al. Lumbar laminectomy alone or with instrumented or noninstrumented arthrodesis in degenerative lumbar spinal stenosis: Patient selection, costs, and surgical outcomes. *Spine* 1997; 22: 1123–1131
13. Lane WA. Case of spondylolisthesis associated with progressive paraplegia: laminectomy. *Lancet* 1893; I: 991
14. Lee CK, Rauschnig W, Glenn W. Lateral lumbar spinal canal stenosis: Classification, pathologic anatomy and surgical decompression. *Spine* 1980; 13: 313–320
15. Lemaire JJ, Sautreaux JL, Chabannes J et al. Lumbar canal stenosis: retrospective study of 158 operated cases. *Neurochirurgie* 1995; 41: 89–97
16. Miller JAA, Schmatz C, Schultz AB. Lumbar disc degeneration: correlation with age, sex and spine levels in 600 autopsy specimens. *Spine* 1988; 13: 173–178
17. Ooi Y, Mita F, Satoh Y. Myeloscopic study on lumbar spinal canal stenosis with special reference to intermittent claudication. *Spine* 1990; 15: 544–549
18. Porter RW. Vascular compression theory. In: Gunzburg R, Szpalski M (eds) *lumbar spinal stenosis*. Philadelphia: Lippincott Williams & Wilkins, 2000: 159–162
19. Sato K, Kikuchi S. Clinical analysis of two-level compression of the cauda equina and the nerve roots in lumbar spinal canal stenosis. *Spine* 1997; 22: 1898–1903

INDIVIDUELLES DESIGN / CUSTOMIZED DESIGN

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You are you looking for an customized design?

Our Endo-Kerrison handles are from a certain volume available in various coloured silicone coatings! Just contact us to create your customized design offer! Your team at ORTHOMEDICAL GmbH