



**NEW K-WIRE DISTRACTORS**

## K-WIRE DISTRACTOR, SINGLE DRILL GUIDE

---

**CAUTION! FOR USE MAXIMUM K-WIRE DIAMETER = DIA. 1.6MM**



11-010.18 K-Wire Distractor 18cm, with Ratchet,  
Max. Ø 1.8mm, single Drill Guide,  
Spreading-Expansion Max. 55mm

## K-WIRE DISTRACTOR, DOUBLE DRILL GUIDE

**CAUTION! FOR USE MAXIMUM K-WIRE DIAMETER = DIA. 1.6/2.0MM**



11-011.18 K-Wire Distractor 17.5cm, with Ratchet,  
Max. Ø 1.8mm/2.2mm, double Drill Guide,  
Spreading-Expansion Max. 65mm

